

Quantum is an elite compact laminate that is part of Novalab's benchtop material range.

"An affordable, high quality benchtop solution."



Properties

- Has an extremely durable and chemical resistant surface
- Completely Impervious to moisture ingress.
- No applied edge lippings.
- Supplied with capillary groove to underside if required.
- See Quantum Properties sheet for chemical and stain resistance information.
- Sheet sizes are 3600 x 1500 x 16mm thk.
- Stock colours held here in New Zealand.

Applications

- Designed specifically for use in standard laboratories.
- Perfect for benchtops and splashbacks in high functioning work and lab areas.
- Can be installed above steel framework such as our System 2000, or above any joinery solution - See BioPro or Novatek.



Environmentally Sustainable



Complete Water Resistance



High Impact Resistance



Cost Effective



10 Year Warranty



Resistant to Chemicals

View the colour range at: www.novalab.co.nz

