

Novalab is New Zealand's only registered distributor for Trespa Toplab Plus
The WORLD'S leading compact laminate for benchtops in laboratories.

"The ultimate educational benchtop material"



Properties

- Trespa Toplab plus is deemed to be the perfect solution for high performance, long lasting science benchtops.
- High level of scratch and graffiti resistance.
- Impervious to moisture ingress.
- The world's only lab benchtop made with chopped strand technology, giving unrivalled core resistance to delamination.
- No applied edge lippings.
- Sheet sizes are 3000 x 1500 & 2550 x 1800 x 16mm thk.
- Stock colours held here in New Zealand.

Applications

- Perfect for all educational laboratory environments.
- Perfect for benchtops and splashbacks in technicians rooms and prep rooms.
- Can be installed above steel framework such as our System 2000, or above any joinery solution - See BioPro or Novatek.



Environmentally
Sustainable



Complete Water
Resistance



High Impact
Resistance



Cost
Effective



10 Year
Warranty



Resistant to
Harsh Chemicals

View the colour range at: www.novalab.co.nz

