

BioPro is Novalab's compact laminate joinery range, specifically designed for use in laboratory areas.

"The perfect solution for durable and water resistant joinery."



Properties

- Joinery made from compact laminate which will not support microbiological growth.
- Impervious to moisture ingress.
- No applied lippings.
- See our BioPro Properties sheet for chemical and stain resistance information.
- Large colour range in stock here in New Zealand.
- Meets the requirements of EN4382 & NEMA LD3-2000

Applications

- Designed specifically for use in PC2 and PC3 laboratories
- Select from our standard joinery range or talk to us about your bespoke requirements.
- Perfect joinery for wet benches, chem labs, micro labs, store areas, splashbacks, safety stations and shelving.
- Designed to work with Novalab's System 1000 / 2000 2032 steel support framework.



Environmentally Sustainable



Complete Water Resistance



High Impact Resistance



Cost Effective



5 Year Warranty



Resistant to Chemicals



View the colour range at: www.novalab.co.nz